

# 941 SENSORS AND ACTUATORS

---

## 1ST & 2ND PRESENCE BUTTON

Code : **9410037**

Model : **LLC-131**

### Description

Room peripheral with two buttons: a presence button with a green LED for the nurse, and a presence button with yellow light for their assistant. This peripheral allows the nurse to indicate the start and/or end of their presence in the room by activating the green present/away button, while the green LED allows to identify the presence status of the nurse in the room. Likewise, their assistant indicates their presence using the yellow present/away button. Wall or surface mount.



9410037

LLC-131

**Especificaciones**

Botón de presencia de enfermera	Color botón: verde (RAL 6018, verde amarillento)
	Color indicador led: verde
	Icono: persona 1
Botón de presencia de auxiliar	Color icono: blanco (RAL 9003, blanco señales)
	Color botón: amarillo (RAL 1007, amarillo narciso)
	Color indicador led: amarillo
Botón de presencia de auxiliar	Icono: persona 12
	Color icono: blanco (RAL 9003, blanco señales)
	Color botón: blanco (RAL 9003, blanco señales)
Alimentación	24 V
Consumo (reposo)	10 mA
Consumo (máximo)	40 mA
Dimensiones (anchura x altura x profundidad)	80x 80x 30 mm
Peso	80g
Material	ABS / Policarbonato
Color (embellecedor)	Blanco (RAL 9001, blanco crema)
Temperatura de funcionamiento	5°C a 55°C
Protección	IP 42
Montaje	Empotrado (con ACC-100 / ACC-101) Superficial (con ACC-110)
Normativa (Seguridad)	EN 60950-1:2006
Normativa (Emisiones electromagnéticas)	EN 61000-3-3