

941 SENSORS AND ACTUATORS

BED UNIT WITH 2 BUTTONS & CONNECTOR

Code : **9410258**

Model : **LLC-409**

Description

Room peripheral that has a DIN8 connector for the connection of the handset (LLC-71x), a call button with a red LED and a nurse presence button.

This peripheral allows the patient to place call with either the pear-type button or the red call button. Once initiated, the red LED indicator allows an easy identification of the bed that originated the call and of its status. The green button is used to either indicate that the nurse is present, or to cancel a call. Wall mount.

