

941 SENSORS AND ACTUATORS


BATHROOM CALL UNIT WITH CANCELLATION

Code : **9410260**

Model : **LLC-421**

Description

Bathroom peripheral with a 2m pull-cord to place calls, a call button with a red LED and a cancellation button. It allows the patient to launch a call by pulling the pull-cord or pushing the red call button. Once the patient call is initiated, the red LED indicator allows an easy identification of the status of the call. This product has an anti-strangle mechanism so that the pull-cord can be split into 2 parts, avoiding accidents. The pull-cord has a red ball at its lower end that allows the patient to easily locate it. Wall mount.



9410260

LLC-421

Specifications

Call pull-cord	Length: 2m red plastic
	Anti-strangulation mechanism 6Kg
	Pull-cord clip color red (RAL 3020, red traffic)
Nurse Call Button	Button color: red (RAL 3020, red traffic)
	LED indicator color: red
	Icon: Nurse Color icon: white (RAL 9003, white signal)
Cancel Button	Button color: (RAL 3020, red traffic)
	Icon: Cancellation
	Color icon: white (RAL 9003, white signal)
Power supply	24 V
Consumption (idle)	2,5 mA
Consumption (maximum)	10 mA
Dimensions	75x120x15 mm
Weight	75g
Housing	ABS + Polycarbonate
Colour	White (RAL 9003, white signal)
Operating temperature	5..55 °C
Protection index	IP 42
Mounting	Recessed Mounting (rectangular box with 4 screw- not available in ALCAD) Superficial (rectangular box with 4 screw- not available in ALCAD)
Safety conformity	EN 60065:2014
EMC conformity	EN 61000-6-3:2007/ A1:2012, EN 61000-6-1:2007